

The removable denture procedure developed in a digital workflow.

## Optima Protocol: Partial

### OPTIMA

REMOVABLE PROSTHETICS PROCEDURE IN FEWER VISITS

Optima is the removable denture solution developed in a digital workflow. In just a few visits, you will get detailed information about your patient's mouth and the work needed, as well as determine which design suits them best, ensuring a natural and beautiful finish.

### BENEFITS



Digital



Accuracy and quality



Fewer visits





Aesthetics

### PACKS

1 — BASIC PATIENT OPTIMA PACK

Contents:

- Optima permanent prosthesis.
- Medicinal mouthwash for better oral protection.
- Special brush for removable prostheses.



In accordance with our commitment as a socially responsible company, we've created sustainable packaging by minimising materials and using recycled and sustainable elements.



# Optima Protocol: Partial



	1ST VISIT	2ND VISIT (WITHOUT OCCLUSION)	3RD VISIT
PROTOCOL ELEMENTS	ANATOMICAL IMPRESSION	OCCLUSION RECORD	FINISH
PROTOCOL ELEMENTS		OPTIMA ELEMENTS – Joint model.  CLINIC MATERIALS – Occlufast silicone.	OPTIMA – Upper or lower Optima partial prosthesis.
PROTOCOL STEPS	1a→ Mouth scan with validated iOS. 1b→ Mouth impression (alginate). 2→ Create prescription* in Corus Link. 3a→ Upload STL files of anatomical impression to Corus Link. 3b→ Send anatomical impression to your Corus Lab. 4→ Upload images to Corus Link.	1→ Check the joint model. 2→ Record the occlusion. 3a→ Use a validated iOS device to scan the occlusion in the patient's mouth with the remaining teeth visible. 3b→ Check the occlusion record in the models. 4→ Update prescription in Corus Link. 5→ Upload images to Corus Link. 6a→ Upload STL files of anatomical impression to Corus Link. 6b→ Send the occlusion record to your Corus Lab.	1→ Validate fit and retention. 2→ Check aesthetic and functional aspects. 3→ Upload images of final result to Corus Link.
ANALOGUE			
DIGITAL			
REQUIREMENTS & MATERIALS	NOTE - If the occlusion with the remaining teeth is not visible, a joint model will be required to correctly take the occlusion records. 2nd VISIT  <input type="radio"/> Send STL. <input type="radio"/> Send impression. <input type="radio"/> Upload prescription.	<input type="radio"/> Send STL. <input type="radio"/> Send the occlusion record. <input type="radio"/> Upload prescription.	<input type="radio"/> Upload final videos and images.
	iOS: Trios 3, 4; Itero; Primescan; Medit 700i; Carestream; Vivascan; Heron. Materials: Alginate.	iOS: Trios 3, 4; Itero; Primescan; Medit 700i; Carestream; Vivascan; Heron.	