CHECKLIST – GUIDED SURGERY

(What do I need on the day of the surgery?)

Language:

ENGLISH

ESPAÑOL

DUTCH

FRANCAIS

INSTANT BOX

	ITEM	DESCRIPTION / REQUIREMENT
☐	Stackable Metal Guide & Pin Guide	The Stackable Metal Guide is used to accurately position all subsequent guides within the guided surgery workflow. It is initially positioned using the Pin Guide and secured to the bone with fixation pins.
☐	Surgical Guide	The Surgical Guide is the one used to carry out the entire drilling protocol and implant placement protocol according to the approved treatment plan.
☐	Provisional Prosthesis & Positioning Guide	The Positioning Guide is used to accurately seat the Instant provisional restoration and capture the Instant abutments.
☐	Abutments	Includes all components required for provisional prosthesis placement: Instant resin abutments, click and clinical screws, and transepithelial abutments.
☐	Surgical Documentation	Includes the Guided Surgery Report, Drill Report (if applicable), and printed or digital access to the treatment planning files.

REQUIRED MATERIAL IN CLINIC

	ITEM	DESCRIPTION / REQUIREMENT
☐	Surgical Kit	The Guided Surgery Kit must correspond to the approved one agreed upon during the initial planning phase.
☐	Resin	Auto-polymerizing resin is used to capture and secure the provisional prosthesis to the Instant abutment cylinders. Recommended: Bredent QU Resin.
☐	Fixation Pins (often included in the surgical kit)	Fixation pins are used to secure the Stackable Metal Guide to the bone, ensuring guide stability throughout the surgical procedure.
☐	Drills for Fixation Pins (often included in the surgical kit)	Dedicated drills used to prepare the osteotomies for fixation pin placement once the Stackable Metal Guide has been positioned.
☐	Screwdriver tip for transepithelial abutments	• IPD/DAS/Medentika: Hex 1.25 • Straumann: Torx • Archimedes: Klockner Star screwdriver
☐	Toolbox	Required for provisional restoration capture and finishing procedures.
☐	PTFE Tape (Teflon)	Used to seal the abutment access channels during provisionalization, preventing resin ingress during the provisional capture procedure.
☐	Dental Dam	Used during provisionalization to control excess resin and prevent apical resin flow. It maintains a clean operative field and facilitates accurate seating of the provisional restoration, verified by a tactile “click”.

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