

Our commitment is to innovate and offer the most advanced dental health solutions that put digital technology at the service of wellbeing.

Avant-garde on human scale.

Corus Product

# Right

Small movements,  
smiles in place.

RIGHTDIPEN\_REV20230222



#### CHANGE THE SWITCH

Right is an aligner system specially designed for prosthodontists. With slight movements, they can correct malpositions and achieve the best possible cosmetic treatments.

#### PLAN AT YOUR SERVICE

The treatment has been designed so that the clinician can carry it out without complications, without having to refer or unnecessarily lengthen the process, with the support of the Right Planning Centre, made up of expert orthodontists.

#### SIMPLE, CLEAR, AGILE

With digital planning, protocols are clear, easy to understand and easy to operate. Because everything is planned in advance, there are no complications. All you have to do is follow the instructions.



#### INVISIBLE ALIGNERS

Right is the first cosmetic orthodontic treatment that is 100% controlled by your dentist. Consisting of invisible aligners that are digitally made to measure, it corrects your smile without affecting your day-to-day life.

#### HOW IT WORKS

You wear each aligner for 10 days before replacing it with the next one. Ideally, you should wear it for around 22 hours a day (no less than 20 hours daily) to get the best results. It works by exerting pressure on specific points of your teeth, causing them to gradually move towards the desired position.

Depending on the complexity, the duration of the treatment can be from one month (Rightmonth) to a maximum of four trimesters (Rightrimesters). At the end of each trimester, there is a digital check-up so that the dentist can ensure the treatment is progressing as planned. Everything under control, at all times.

MORE INFORMATION  
Ask your dentist.



#### IT MOVES WITH YOU

We give you your first aligners along with a storage case, a special maintenance kit and a travel pouch that lets you take Right with you at all times.

#### Your smile, in place

Through small, carefully calculated movements, your teeth are gradually straightened and you get the smile you've always wanted. Because beauty is in the details.

#### Total discretion

Wearing invisible aligners means that you're the only one who knows you're wearing them and you can go about your life as normal. The only thing you have to do is take them out during meals. Both imperceptible and painless.

#### For a full smile

Aligning your teeth with Right is the perfect basis for future whitening or veneer treatments that will bring your smile to perfection.

#### In confidence

The whole treatment is in the hands of your dentist. Because you can trust them and they know your smile best, everything is easier and simpler and you spend less time in the chair.

RRight is a custom-made medical device, manufactured in accordance with the requirements of the European Directive 93/42/CEE in health centres authorised by the competent authorities.